



7105 Northland Terrace
 Brooklyn Park, MN 55428
 Ph# 1-888-538-6872
 Fax# 763-582-0442
 Amy Jorgensen-Inside Sales
AmyJ@LevoUSA.com
www.levousa.com

USA Quote / Order Form LCEV Kid / Jr



Dealer Information

Contact _____ Fax# _____
 Company _____ Ph# _____
 Address _____ Email _____
 City _____ Date _____
 State/Zip _____

Client Information

Name _____
 Height: _____ Weight: _____
 Accurate weight is important for selection
 of proper power components

Please answer the following questions accurately and completely. Identifying the exact measurements within the diagrams allows for correct configuration and final adjustments for the client. Please do not provide measurements from client's existing chair or cushion! Please be

Upper Leg Length- Seat Depth



While sitting, measure the distance between the rear point of the bottom and the hollow of the knee

11" XS - S
 ←————→ 17.75"
 9.5"
 10"
 11.5"
 12"
 13"
 13.75"

Hip Size- Seat Sling Width



While sitting, measure the outermost points of the hips

11.75" ←————→ 13.75"
 11.75"
 13.75"


Lower Leg Length



While sitting, measure from the leg to bottom of heels with shoes on. Do not use a cushion






7.5" ←————→ 16"
Footrest position
 7.5"
 8.25"
 9"
 9.5"
 10.5"
 11.25"
 11.75"
 12.25"
 13.25"
 14"
 14.5"
 15.25"
 15.75"








Backrest Height






Measure the distance between the seating surface and the lower shoulder blade bone

"R" _____

Armrest Height		Bitte
	Distance between seat surface and bottom of forearm	7" ← _____ → 9.75"
Headrest Height Bitte		
	From the top of seat surface to mid-point back of head	"H" _____
LCEV Kid/JR Base- K0005		
Manual wheelchair with following features as standard: auto stability lock in all seating and standing positions; adjustable frame depth; armrest round with built-in control switch; fold down and angle adjustable backrest; side guards; extra strong knee support; basic frame color; quick release front and rear wheels; rear air filled tires with aluminum handrims; standard locking brakes and solid rubber front casters Crash tested per ISO 7176/19 ANSI/RESNA WC/Vol. 1-Section 19		1,000.00 <input type="checkbox"/>
Powered Standing Function- E2301 Must be ordered		11,875.00 <input type="checkbox"/>
Hand Held Control <input type="checkbox"/> Switch in Arm- Right <input type="checkbox"/> Switch in Arm-Left <input type="checkbox"/>		
Operation		
Right <input type="checkbox"/> Left <input type="checkbox"/>		
Frame Colors -- Check preference		Special Colors -- Additional charge \$490.00 <input type="checkbox"/>
Standard Colors -- No charge		
Traffic Red _____	Red Metallic _____	Orange Reflective _____
Zinc Yellow _____	Green Metallic _____	Lemon Reflective _____
Marine Blue _____	Blue Metallic _____	Silver Reflective _____
Traffic Black _____	Silver Metallic _____	
Platinum _____	Polished Chrome _____	\$950.00 <input type="checkbox"/>
Chest Support Belt- E0960		\$120.00 <input type="checkbox"/>
Must be ordered		
Armrest / Control Options		
Standard- Armrest Round with Control Switch- K0108		No charge <input type="checkbox"/>
Option- Armrest round with hand control		No charge <input type="checkbox"/>
Folding Down Stroller Bar - Width Adjustable		\$823.00 <input type="checkbox"/>
	Option- Up / Down Hand Held Device for Attendant Price includes Stroller Bar and Hand Held Device	\$1,218.00 <input type="checkbox"/>
Dual Switch for Table Control- K0108		\$1,160.00 <input type="checkbox"/>
	Stand-up operation with push buttons for table Creates a jack input for the dual switch Table tray not included	
Control Switch Extension- K0108		Vertical for round armrest \$102.00 <input type="checkbox"/>
	Horizontal \$30.00 <input type="checkbox"/>	

Backrest Options			
Choose Backrest Angle			
Angle Adjustable Backrest Set in Position 1- 85°		No charge	<input type="checkbox"/>
Angle Adjustable Backrest Set in Position 2- 90°		No charge	<input type="checkbox"/>
Angle Adjustable Backrest Set in Position 3- 95°		No charge	<input type="checkbox"/>
Omit Backrest		No charge	<input type="checkbox"/>
Client can choose alternative backrest- LCEV comes standard with dual cane backrest supports- Please specify preferred back cane height			
Tension Adjust Sling Backrest- K0108		\$135.00	<input type="checkbox"/>
Choose backrest cane height	12"	No charge	<input type="checkbox"/>
	13"	No charge	<input type="checkbox"/>
	13.75"	No charge	<input type="checkbox"/>
	14.5"	No charge	<input type="checkbox"/>
	15.25"	No charge	<input type="checkbox"/>
	16"	No charge	<input type="checkbox"/>
Backrest Extension 6"- K0108		\$335.00	<input type="checkbox"/>
	Used only on Tension Adjust Sling Backrest that are higher than 16"		
V-Trak Standard & Axxis Backrest- E2620			
	V-Trak Standard Backrest -- 11" High - Can only be used on 30cm wide chair	\$837.00	<input type="checkbox"/>
	V-Trak Standard Backrest -- 12" High	\$837.00	<input type="checkbox"/>
	V-Trak Axxis Backrest -- 12" High	\$1,012.00	<input type="checkbox"/>
	V-Trak Axxis Backrest -- 16" High	\$1,113.00	<input type="checkbox"/>
V-Trak Backrest Adjustable Hardware- E1028		\$292.00	<input type="checkbox"/>
Must be ordered with V-Trak Backrest			
V-Trak Thoracic Accessories			
	V-Trak Thoracic Lateral Pads Large -- 5" x 7"-- Right & Left- E0956	\$272.00	<input type="checkbox"/>
	V-Trak Thoracic Lateral Pads Medium -- 4.25" x 6"-- Right & Left- E0956	\$272.00	<input type="checkbox"/>
	V-Trak Thoracic Adjustable Hardware- E1028	\$426.00	<input type="checkbox"/>
Must be ordered with V-Trak Lateral Pads			
V-Trak Headrest- E0955			
	V-Trak Headrest Winged	\$312.00	<input type="checkbox"/>
	V-Trak Headrest Simple	\$287.00	<input type="checkbox"/>
	V-Trak Headrest Adjustable Hardware- E1028	\$231.00	<input type="checkbox"/>
Must be ordered with V-Trak Headrest			
Rear Wheel Options & Accessories			
	Rear wheel 16", castor 3"	Kid 33 (11.5")	No charge <input type="checkbox"/>
	Rear wheel 20", castor 5"	Junior 40 (14.5")	No charge <input type="checkbox"/>
	Rear wheel 22", castor 5"	Junior 44 (16")	No charge <input type="checkbox"/>
Corresponding seat to floor height without cushion			
Rear Wheels Spinergy- K0108			
	Light weight, improved driving performance, easy to clean Includes mounted aluminum handrims Can only be used on 22" rear wheels Rim color is black. Spoke colors available: black; red & white	18 Spoke (Pair)	\$697.00 <input type="checkbox"/>
Brake Lever Extension- 2.4"- E0961			
	Simpler accessibility and less strength necessary Standard pneumatic or solid rubber rear wheels only	Right Left	\$75.00 <input type="checkbox"/> \$75.00 <input type="checkbox"/>

Knee Support Options			
Rectangular- K0108		Conical- K0108	
	2" x 4" \$35.00 <input type="checkbox"/> 3" x 4" \$35.00 <input type="checkbox"/> 3" x 6" \$35.00 <input type="checkbox"/>		2" x 4" x 4" \$44.00 <input type="checkbox"/> 3.5" x 5.5" x 4" \$44.00 <input type="checkbox"/> 5" x 7" x 4" \$44.00 <input type="checkbox"/>
Knee Support Height Adapter- K0108			
Ensures proper alignment with knees for clients with longer lower legs		Pair	\$298.00 <input type="checkbox"/>
Footplate Options			
Choose Footplate Angle			
Angle Adjustable Footplate Set in Position 1- 80°			No charge <input type="checkbox"/>
Angle Adjustable Footplate Set in Position 2- 90°			No charge <input type="checkbox"/>
Angle Adjustable Footplate Set in Position 3- 100°			No charge <input type="checkbox"/>
1-Piece Fixed -- No Flip-Up to Rear- K0040			\$140.00 <input type="checkbox"/>
			
Footrest Accessories			
Seat Depth / Footrest Quick Adjustment- K0108			\$118.00 <input type="checkbox"/>
	Allows for quick adjustment of seat depth and footrest height		
Calf Strap- K0038			\$73.00 <input type="checkbox"/>
Must be ordered			
Foot Guides- K0108			<input type="checkbox"/>
	Allows individual foot adjustment to accommodate client needs while increasing stability	Small Right	\$213.00 <input type="checkbox"/>
		Small Left	\$213.00 <input type="checkbox"/>
		Large Right	\$233.00 <input type="checkbox"/>
		Large Left	\$233.00 <input type="checkbox"/>
Non-Slip Footrest Surface- K0108			\$129.00 <input type="checkbox"/>
Accessories			
Foam Seat Cushion with Seat Cushion Cover- E2601			
Additional Seat Cushion Cover- E2619			
	Seat Cushion		\$397.00 <input type="checkbox"/>
	Additional Seat Cushion Cover		\$219.00 <input type="checkbox"/>
Backpack			\$124.00 <input type="checkbox"/>
	Yellow <input type="checkbox"/>	Red <input type="checkbox"/>	
	Black <input type="checkbox"/>	Blue <input type="checkbox"/>	
Maxgrip Handrims- K0108			
	Only for 22" rear wheels		
	Handrims- rubber covered		
	Original	\$578.00 <input type="checkbox"/>	
	Ergo Para	\$915.00 <input type="checkbox"/>	
Table Tray / Fixed- E0950			\$282.00 <input type="checkbox"/>
	Supporting hardware installed on the armrest. Easily attached and removed Control Switch Extension or Hand Held Device needed to run chair function with Table Tray attached		

Anti-Tip for Kid- E0971



Right	\$103.00	<input type="checkbox"/>
Left	\$103.00	<input type="checkbox"/>

Anti-Tip for Junior- E0971



Right	\$106.00	<input type="checkbox"/>
Left	\$106.00	<input type="checkbox"/>

Spoke Guard-Worms- K0065

\$186.00



Plexiglass shields to protect hands and fingers

Chair Order Policy Standardausführung der Training-Dive-Option. Ermöglicht

Contact LEVO USA at 1-888-LEVOUSA with any questions regarding how to accurately complete this form.

LEVO USA is not responsible for configuration or size discrepancies resulting from inaccurate or incomplete information provided on the order form.

Order confirmations will be sent to you. Please review and sign your agreement/confirmation authorizing submission of the order as stated. Return the signed copy to LEVO USA at fax number 763-582-0442 or via email at: amyj@danetechnologies.com.

All chairs are built to specifications provided, therefore, all sales are final. No returns or exchanges. Your signature acknowledges your understanding and acceptance of this no return policy. Customer is responsible for full payment of chair shipped based on specifications provided on confirmed order form.

The HCPC codes identified herein are our best interpretation of the match between chair, accessories and the appropriate code. This coding interpretation has not been verified by DME PDAC nor should our interpretation replace the customer's own responsibility for submitting accurate information to third party sources.

Please note: All prices and specifications provided on form are subject to change without prior notice. All prices for options or accessories on order form are only valid if ordered with the chair or prior to completed production of unit. Any chair modifications, options or accessories required for future use or maintenance of chair will carry different pricing and warranty terms.

Warranty Policy

Your LEVO product is guaranteed from the date of purchase for:

Two years covering all material and manufacturing defects of mechanical parts.

LEVO USA will not repair or replace free of charge any part or parts found to be defective due to abuse, misuse or lack of maintenance.

