



Dealer: _____ Phone: _____ Client: _____

Contact: _____ Email: _____

Address: _____ PO#: _____ Date: _____

City/State/Zip: _____

LEVO USA
 7105 Northland Terrace
 Brooklyn Park, MN 55428
 888.538.6872
 888.Levousa
 763.582.0442 FAX
 sales@levousa.com
 AmyJ@levousa.com
 www.levousa.com




Purchasing through (funding source): _____

Standard Base Configuration \$9,595

- Manually Propelled Wheelchair with Manual Standing System
- Maximum Client Weight: 265 lbs.
- Gas Strut Assist to Stand
- Lateral Lifting Arms/Sideguards
- Tension Adjustable Seat & Back Upholstery
- Fold Down, Angle Adjustable Backrest
- Rear Wheel - 24"x1" high pressure pneumatic, hub black, spoke silver
- Front caster - 5"x1" solid rubber, composite wheel
- Standard colors: platinum, red, yellow, blue or black
- One Piece Footplate - Fixed
- 0 Degree Rear Wheel Camber

Anatomical Dimensions

Weight:	
Height:	
Hip Width: ①	
Upper Leg Length: ②	
Lower Leg Length: ③	
Seat to Elbow: ④	
Chest Width: ⑤	
Chest Depth: ⑥	
Floor to top of knee: ⑦	
Seat to Top of Shoulder: ⑧	



① Hip width > Seat width

Hip width:

Seat width:	
13.75"	Option
15"	Option
15.75"	Option
16.5"	Option
17.75"	Option
19"	Option

\$ 1,035

② Upper leg length/Seat depth (upholstery)

Upper leg length:

There are two frame lengths for the LAE (S-M & L-XL)
 Each offers a range of seat depths utilizing the adjustable upper and lower seat rails.

When fit correctly, the *Seat depth (upholstery)* is approximately 2.5" less than the *Upper leg length* measurement.

Seat depth (upholstery):		
S-M Frame	11.5"	Option
	12"	Option
	13"	Option
	14"	Option
	15"	Option
L-XL Frame	16"	Option
	17"	Option
	18"	Option
	18.5"	WAMB \$ 199
	19.5"	WAMB \$ 199
	20"	WAMB \$ 199
21"	WAMB \$ 199	

③ Lower leg length/Seat to footplate

The *Lower leg length* determines the *Seat to footplate* setting, possible wheel/caster combinations and front seat heights.

Lower leg length:

Please take the cushion into consideration when choosing the *Seat to footplate* setting.

Seat Height Front: ▶	19"	20"	21"	22.5"	Rear seat height is approximately 1.5" lower than front
Seat to footplate: ▼					
	22" rear wheel 5" front caste	24" rear wheel 5" front caste	24" rear wheel 5" front caste	26" rear wheel* 6" front caster*	
14"	<input type="checkbox"/>				
15"	<input type="checkbox"/>				
15.75"	<input type="checkbox"/>	<input type="checkbox"/>			
16.5"	<input type="checkbox"/>	<input type="checkbox"/>			
17.25"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
18"		<input type="checkbox"/>	<input type="checkbox"/>		
19"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
19.75"			<input type="checkbox"/>	<input type="checkbox"/>	
20.5"			<input type="checkbox"/>	<input type="checkbox"/>	
21.25"			<input type="checkbox"/>	<input type="checkbox"/>	
22"			<input type="checkbox"/>	<input type="checkbox"/>	

*see Wheel options on page 4

Client weight

Client weight:

Gas strut selection chart

Depth: ▶	11.5"	12"	13"	14"	15"	16"	17"	18"	18.5"	19.5"	20"
Weight: ▼											
88 lbs. to 132 lbs.	<input type="checkbox"/> 400n	<input type="checkbox"/> 400n	<input type="checkbox"/> 400n	<input type="checkbox"/> 400n	<input type="checkbox"/> 400n	<input type="checkbox"/> 400n	<input type="checkbox"/> 400n	<input type="checkbox"/> 400n	<input type="checkbox"/> 400n	<input type="checkbox"/> 400n	
132 lbs. to 176 lbs.			<input type="checkbox"/> 400n	<input type="checkbox"/> 500n	<input type="checkbox"/> 500n	<input type="checkbox"/> 500n	<input type="checkbox"/> 500n	<input type="checkbox"/> 500n	<input type="checkbox"/> 500n	<input type="checkbox"/> 500n	
176 lbs. to 220 lbs.				<input type="checkbox"/> 500n	<input type="checkbox"/> 600n	<input type="checkbox"/> 600n	<input type="checkbox"/> 600n	<input type="checkbox"/> 600n	<input type="checkbox"/> 600n	<input type="checkbox"/> 600n	<input type="checkbox"/> 600n
220 lbs. to 264 lbs.							<input type="checkbox"/> 600n	<input type="checkbox"/> 700n	<input type="checkbox"/> 700n	<input type="checkbox"/> 700n	<input type="checkbox"/> 700n

Frame color

Black	Std.		•
Blue	Option		•
Platinum	Option		•
Red	Option		•
Yellow	Option		•
Blue metallic	Option		\$ 595
Green metallic	Option		\$ 595
Red metallic	Option		\$ 595
Silver metallic	Option		\$ 595
Orange reflex	Option		\$ 695
Lemon reflex	Option		\$ 695
Silver reflex	Option		\$ 695
Aluminium polished	Option		\$ 995
Traffic white	Option		\$ 895
Other RAL color	Option		\$ 895

RAL: 9016

Standing Release

Bilateral	Std.		•
Right	Option		•
Left	Option		•

Backrest


Backrest - canes with tension adjustable sling and cushion

	Upholstery height:		
12.25"	WAMD		\$ 111
13.00"	WAMD		\$ 111
13.75"	WAMD		\$ 111
14.5"	WAMD		\$ 111
15.5"	WAMD		\$ 111
16"	WAMD		\$ 111

Backrest - V-Trak "Dynamic Stability" with internal adjustable straps and adjustable laterals (Flexes to absorb energy, .25" padded/breathable upholstery (must order mounting hardware)

V-TRAK Axxis with parallel sides

Parallel Sides with 2-3 inches of contour depth



Internal adjustable spine profiling straps





Internal adjustable lateral support

Width ▶	M	L	XL
Adjustable:	Inside Width 9.5"-12"	Inside Width 11"-13.5"	Inside Width 13.5"-15.75"
① Height 12"	\$695 <input type="checkbox"/> E2620	\$695 <input type="checkbox"/> E2620	\$750 <input type="checkbox"/> E2620
② Height 15.75"	\$695 <input type="checkbox"/> E2620	\$695 <input type="checkbox"/> E2620	\$750 <input type="checkbox"/> E2620
③ Height 21.5"	\$750 <input type="checkbox"/> E2620	\$775 <input type="checkbox"/> E2620	\$795 <input type="checkbox"/> E2620



V-TRAK Axxis Deep Contour (Lateral crest is always 22cm from bottom)

Deep contour with 5-7 Inches of contour depth




Internal adjustable spine profiling straps

Internal adjustable lateral support - Lateral adjustment for front of backrest

Width ▶	M	L	XL	XXL
Adjustable:	Inside Width 9.5"-12"	Inside Width 11"-13.5"	Inside Width 13.5"-15.75"	Inside Width 14"-21.5"
④ Height 21.5"	\$775 <input type="checkbox"/> E2620	\$795 <input type="checkbox"/> E2620	\$825 <input type="checkbox"/> E2620	\$850 <input type="checkbox"/> E2620
⑤ Height 23.5"	\$775 <input type="checkbox"/> E2620	\$795 <input type="checkbox"/> E2620	\$825 <input type="checkbox"/> E2620	\$850 <input type="checkbox"/> E2620



Mounting Hardware - Adjustable



1	Lite - tube style			\$	292
	Additional set (recommended for standing)			\$	292
2	Original - flat style (recommended for chair back widths over 19")			\$	292
	Additional set (recommended for standing)			\$	292

Backrest options




3	Standard Backrest - 6" tube extension for sling back (standard back height is 16", this item extends it to 22")			\$	369
	V-Trak backrest - 8" extension			\$	543

Backrest angle

85°	Option		•
90°	Std.		•
95°	Option		•

Seat Cushion - memory foam with cover



Seat width	13.75"	15"	15.75"	16.5"	17.75"	19"			\$
Seat depth									
Matches seat depth unless otherwise specified	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			397

Additional Seat Cushion Cover

Standard	SESND1050		\$	219
Incontinent	SESND1055		\$	219

Supports



1	Knee support - Padded and reinforced			\$	172
2	Rectangular knee wedge	2" x 4"	SEFND0035	\$	35
		2.5" x 4"	SEFND0035	\$	35
		3" x 4"	SEFND0035	\$	35
		2" x 6"	SEFND0035	\$	35
		3" x 6"	SEFND0035	\$	35
3	Conical knee wedge	2"- 4" x 4"	SEFND0035	\$	35
		3.5"- 4.5" x 4"	SEFND0035	\$	35
		5"- 7" x 4"	SEFND0035	\$	35
4	Knee support - Sport			\$	295
5	Knee support - Rigid tube, swingaway	Size S (4" - 4.75")		\$	545
		Size M (5.5" - 6.5")		\$	545
	Calf strap (see pict #1 below)	SEFND0010	X	\$	73
	Chest support belt	SEFND0030	X	\$	120
6	Additional Chest support - BodyPoint Universal Neoprene- 4.5" wide - SP110L			\$	195
7	Shoe Holder	3.75" x 7.25" Small	Right	\$	213
		3.75" x 7.25" Small	Left	\$	213
		4.5" x 8" Large	Right	\$	233
		4.5" x 8" Large	Left	\$	233

Footrest



1	One piece footplate, fixed	SEFND0015		\$	95
2	Two-piece footplate, flip-up			\$	285
3	One piece footplate, flip-up (back) (not available with footrest height adapter)			\$	219
2	Footplate height adapter (available only with one piece footplate, fixed) (footplate will be 1" narrower than normal)	For seat widths 15", 17.75" & 19"		\$	69
		For seat widths 15.75" & 16.5"		\$	225
4	Easy transfer footplate (extended "L" brackets bring footplate closer to the floor)	Footplate Depth 8"		\$	388
		Footplate Depth 10.25"		\$	388
	Detachable legrest/downposts			\$	149
	Non-slip footplate surface			\$	129


Footplate angle

80°			\$	99
90°			\$	99
100°			\$	99

Rear wheel position

There are multiple Axle Blocks available on the LAE allowing modifiable fore/aft positions of the rear axle and wheel. The chair comes delivered with the rear Axle Block (and rear wheel position) that provides the most optimum propelling characteristics for the client. For further information, please contact customer service.

Wheel options



Tire - 26"x1" high pressure pneumatic, hub black, spoke silver with 6" front caster	Pair			\$	325											
Tire - 24" Solid Polyurethane	Pair			\$	159											
1 Spiner wheel - 24" w/black Marathon tire		<table border="1"> <tr> <td>Red</td> <td>Black</td> <td>White</td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	Red	Black	White	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	12 Spoke	Pair		\$	925
Red	Black	White														
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>														
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>														
			18 Spoke	Pair		\$	697									
Not available when using hub lock, drum brakes or one arm drive																
2 Camber				3°			\$	274								
Not available when using hub lock, drum brakes or one arm drive				6°			\$	274								
3 Steel hand rim	Pair					\$	132									
4 Maxgrepp - Original handrim with rubber coating	Pair	WAMO				\$	578									
5 Maxgrepp - Ergo Para CP handrim	Pair					\$	995									
6 Rear wheels - hub lock brake	Pair					\$	550									
7 One arm drive	Right	WAML				\$	1,595									
Available only with standard wheels	Left	WAML				\$	1,595									

Accessories



1 Tray table with hardware - removeable		SETND0005		\$	378
2 Push handles	Pair			\$	175
3 Push handles - 2.75" extension	Pair			\$	338
4 Push handles - detachable, height adjustable 1" - 6"	Pair			\$	349
5 Anti Tips - traditional	Pair			\$	152
6 Anti Tips - swingaway with attendant tip aid	Pair			\$	450
XL Adapter					
Necessary with swingaway anti-tips on extra long (XL) frame	Pair			\$	152
7 Tip Aid - attendant	Pair			\$	85
8 Spoke guard - transparent	Pair	WAMP		\$	186
9 Shoe holder	Pair			\$	426
10 Wheel lock extension - 2.5"	Pair			\$	150
11 Seat depth and footrest quick adjustment kit				\$	129
Includes 10 clamps with post and 4 threaded knobs					
12 Frame reinforcement kit				\$	290
Necessary when e-fix or similar system is mounted					
Anti-tips - swingaway required					
Cannot be used with Hub Lock Brakes or Anti Tips - traditional					
Vehicle tie-down kit				\$	988
Fixed attachment points for chair transport					
Includes frame reinforcement kit					
Anti-tips - swingaway required					
Mandatory if client is transported in vehicle. Chair required to be secured with the Q'Straint 5001-T2 or any restraint system certified by Q'Straint 5001-T2. Occupant must be secured using the Q'Straint Vehicle 3 point or any restraint system certified by ISO 10542/2 standards.					

Remarks