

LEVO combi Jr

The Mid-Wheel power wheelchair with stand-up function and complete seating systems for junior sizes

LEVO combi Jr



LEVO stands for a better quality of life

LEVO.dt⁺
THE EXPERTS IN STANDING

Even our children can stand like giants



simply attend to...



The powered LEVO combi Jr means secure and active support for your child. The extremely agile power wheelchair optimizes participation into all of your child's opportunities and activities. School – Home – Play. Like the adult version, the LEVO combi Jr has multiple power seating systems including the LEVO designed and industry leading standing capability. (Your child can drive the LEVO combi Jr even while in the standing position.) Easy adjustable growth plates, components and a wide range of accessories and options let the LEVO combi Jr smoothly grow with your child for many years.



LEVO Stands for a better quality of life

Highlights

- User-friendly, simple handling
- Smooth seat/standing positioning movement
- Compact maneuverability for tight spaces (turning circle 43 in)
- Safe and secure stand-up operation elevating to any and all angles between sitting and 85° standing
- Adjustable for upper and lower leg, foot position and back angle
- Max. user weight up to 265 lbs
- Wide range and style of options, accessories and control switches
- Multiple color options
- Easy accessibility for maintenance and service updates
- Certificates: TÜV, CERAH, HI, FDA/510(k), Crash Test



Technical Data

Model	combi		combi/Jr
	S-M	L-XL	XS-S
Seat width (SW)	15/17.3/19.7		11.8/13.8/15
Total width (W)	26		26
Total length (L)	43.3		43.3
Total height (incl. Backrest upholstery) (H)	35-40.9		35-40.9
Seat height (with seat cushion)	22-27.9		22-27.9
Seat upholstery depth	11.4-16.5	15.7-20.9	9.4-12.2
Armrest height (from top of cushion)	6.7-9.8		5.9-13
Back height (with standard backrest upholstery)	16.9		16.9
Foot rest height (from top of cushion)	14.6-22.8		9-18.5
Turning circle	43		43
Speed	3.7/5/6.2 mph		3.7/5 mph
Range	22 miles		22 miles
Max. gradient	12° (25 %)		12° (25 %)
Curb climbing ability	3.1		3.1
Max. chair weight	375 lbs		375 lbs
Max. user weight	265 lbs		265 lbs
Transportation dimension	36.6 x 6 x 26 x 26.4		36.6 x 26 x 26.4

All measurements in inch and pound. Discrepancies in the measurements may arise due to wheelchair modifications.

Standard Configuration

- Backrest upholstery
- Flip-up armrests, height adjustable
- Flip-up footrest
- Knee support and chest strap
- Length adjustment backrest and height adjustment armrest
- Manual seat height adjustment
- Powered stand-up function
- Powered driving in sitting and standing
- Powered lying-standing direct function if option reclining included
- Set of batteries and charger
- PGDT controller (P+/R-net)

Options/Accessories

- Swing away joystick
- Lights and indicators
- V-Trak back system
- V-Trak head and thoracic support
- Powered seat lift 5.9 in
- Powered reclining backrest / legrest 90 - 180°
- Foam seat cushion 2.4 in
- Assorted knee supports and footrest wedges
- Chest role
- Table tray
- Skirt guard
- Vehicle tie-down kit
- Push handle bar
- ... many more on request

Colors

5 standard colors (blue, red, yellow, platinum, black) and 190 other color options

