



LEVO ^{of} 
 THE EXPERTS IN STANDING
 USA Quote / Order Form LAE

7105 Northland Terrace
 Brooklyn Park, MN 55428
 Ph# 1-877-538-6872
 Fax# 763-582-0442
 Amy Jorgensen-Inside Sales
amyj@danetechnologies.com
www.levousa.com



Dealer Information

Contact _____	Fax# _____
Company _____	Ph# _____
Address _____	Email _____
City _____	Date _____
State/Zip _____	

Client Information

Name _____

Height: _____ Weight: _____ (265 Lbs Max)
 Accurate weight is important for selection of proper power components

Please answer the following questions accurately and completely. Identifying the exact measurements within the diagrams allows for correct configuration and final adjustments for the client. Please do not provide measurements from client's existing chair or cushion! Please be

Upper Leg Length- Seat Depth



While sitting, measure the distance between the rear point of the bottom and the hollow of the knee

S - M		L - XL	
14" ←	→ 19"	19.5" ←	→ 25"
14" <input type="checkbox"/>		19.5" <input type="checkbox"/>	
15" <input type="checkbox"/>		20.5" <input type="checkbox"/>	
16" <input type="checkbox"/>		21.75" <input type="checkbox"/>	
17" <input type="checkbox"/>		22.5" <input type="checkbox"/>	
17.75" <input type="checkbox"/>		23" <input type="checkbox"/>	
19" <input type="checkbox"/>		25" <input type="checkbox"/>	

Hip Size- Seat Pan Width



While sitting, measure the outermost points of the hips

13.75" ←	→ 19"	13.75" <input type="checkbox"/>
		15" <input type="checkbox"/>
		15.75" <input type="checkbox"/>
		16.5" <input type="checkbox"/>
		17.75" <input type="checkbox"/>
		19" <input type="checkbox"/>


Lower Leg Length



While sitting, measure from the leg to bottom of heels with shoes on
 Do not use a cushion




Lower Leg Length	Footrest position		
14" ←	→ 22"	f 19" r 17.7"	14" <input type="checkbox"/>
		f 20" r 18.5"	15" <input type="checkbox"/>
		f 21.2" r	15.75" <input type="checkbox"/>
		f 22.5" r 20.5"	16.5" <input type="checkbox"/>
			17.25" <input type="checkbox"/>
			18" <input type="checkbox"/>
			19" <input type="checkbox"/>
			19.5" <input type="checkbox"/>
			20.5" <input type="checkbox"/>
			21" <input type="checkbox"/>
			22" <input type="checkbox"/>









Backrest Height

























Measure the distance between the seating surface and the lower shoulder blade bone


"R" _____

Armrest Height	
	Distance between seat surface and bottom of forearm 7" ← → 9.75"
Headrest Height	
	From the top of seat surface to mid-point back of head "H" _____
LAE Base- E2230	
Manual wheelchair with following features as standard: auto stability lock in all seating and standing positions; manual gas spring assisted standing function; adjustable frame depth; fold down and angle adjustable backrest; side guards; extra strong knee support; basic frame color; quick release front and rear wheels; rear air filled tires with aluminum handrims; standard locking brakes and solid rubber front casters.	
Crash tested per ISO 7176/19 ANSI/RESNA WC/Vol. 1-Section 19	
Operation	
Right <input type="checkbox"/>	Left <input type="checkbox"/>
Bilateral <input type="checkbox"/>	
Frame Colors -- Check preference	
Standard Colors -- No charge	
Traffic Red _____ Zinc Yellow _____ Marine Blue _____ Traffic Black _____ Platinum _____	Special Colors -- Additional charge
	\$490.00 <input type="checkbox"/>
Red Metallic _____ Green Metallic _____ Blue Metallic _____ Silver Metallic _____ Polished Chrome _____	Orange Reflective _____ Lemon Reflective _____ Silver Reflective _____
	\$950.00 <input type="checkbox"/>
Chest Support Belt- E0960	
Must be ordered	
\$120.00 <input type="checkbox"/>	
Non-Standard Seat Frame Measurements	
Hip width measurement 19" - 19.75"	\$950.00 <input type="checkbox"/>
Upper leg measurement 20"- 21"- E2203	\$734.00 <input type="checkbox"/>
Upper leg measurement 22"- 24"- E2204	\$1,247.00 <input type="checkbox"/>
Backrest Options	
Choose Backrest Angle	
Angle Adjustable Backrest Set in Position 1- 85°	No charge <input type="checkbox"/>
Angle Adjustable Backrest Set in Position 2- 90°	No charge <input type="checkbox"/>
Angle Adjustable Backrest Set in Position 3- 95°	No charge <input type="checkbox"/>
Omit Backrest	
No charge <input type="checkbox"/>	
Client can choose alternative backrest- LAE comes standard with dual cane backrest supports	
Tension Adjust Sling Backrest- K0108	
\$135.00 <input type="checkbox"/>	
Choose backrest cane height	12" <input type="checkbox"/> 13" <input type="checkbox"/> 13.75" <input type="checkbox"/> 14.5" <input type="checkbox"/> 15.25" <input type="checkbox"/> 16" <input type="checkbox"/>
Backrest Extension 6"- K0108	
\$335.00 <input type="checkbox"/>	
	Used only on Tension Adjust Sling Backrest that are higher than 16"


V-Trak Standard & Axxis Backrest- E2620			
	V-Trak Standard Backrest -- 12" High		\$837.00 <input type="checkbox"/>
	V-Trak Standard Backrest -- 16" High		\$837.00 <input type="checkbox"/>
	V-Trak Standard Backrest -- 24" High		\$852.00 <input type="checkbox"/>
	V-Trak Standard Backrest -- 31.5" High		\$920.00 <input type="checkbox"/>
	V-Trak Axxis Backrest -- 12" High		\$1,012.00 <input type="checkbox"/>
	V-Trak Axxis Backrest -- 16" High		\$1,113.00 <input type="checkbox"/>
	V-Trak Axxis Backrest -- 24" High		\$1,213.00 <input type="checkbox"/>
V-Trak Backrest Adjustable Hardware- E1028			\$292.00 <input type="checkbox"/>
	Must be ordered with V-Trak Backrest		
V-Trak Thoracic Accessories			
	V-Trak Thoracic Lateral Pads Large -- 5" x 7"-- Right & Left- E0956		\$272.00 <input type="checkbox"/>
	V-Trak Thoracic Lateral Pads Medium -- 4.25" x 6"-- Right & Left- E0956		\$272.00 <input type="checkbox"/>
	V-Trak Thoracic Adjustable Hardware- E1028		\$426.00 <input type="checkbox"/>
Must be ordered with V-Trak Lateral Pads			
V-Trak Headrest- E0955			
	V-Trak Headrest Winged		\$312.00 <input type="checkbox"/>
	V-Trak Headrest Simple		\$287.00 <input type="checkbox"/>
	V-Trak Headrest Adjustable Hardware- E1028		\$231.00 <input type="checkbox"/>
Must be ordered with V-Trak Headrest			
V-Trak Backrest Extension 8"- K0108			\$543.00 <input type="checkbox"/>
Used only on V-Trak Backrest without Push Handle Bar			
Rear Wheel Options & Accessories			
	Rear wheel 22", castor 5"	Corresponding seat to floor height without cushion, front (f) and rear (r) measurements	f 19.0" r 17.7" No charge <input type="checkbox"/>
	Rear wheel 24", castor 5"		f 20.0" r 18.5" No charge <input type="checkbox"/>
	Rear wheel 24", castor 5"		f 21.2" r 19.3" No charge <input type="checkbox"/>
	Rear wheel 26", castor 6"		f 22.5" r 20.5" \$295.00 <input type="checkbox"/>
Rear Wheels Solid Rubber Tires- E2220			\$145.00 <input type="checkbox"/>
Includes mounted aluminum handrims			
Rear Wheels Marathon Black- K0108			\$115.00 <input type="checkbox"/>
Pneumatic tires with tubes and inserts Includes mounted aluminum handrims Recommended when ordering Spinergy Rear Wheels			
Rear Wheels Spinergy- K0108			
	Light weight, improved driving performance, easy to clean	12 Spoke (Pair)	\$925.00 <input type="checkbox"/>
	Includes mounted aluminum handrims Cannot be used with Hub Locking Brakes or E-Motion Wheels Rim color is black. Spoke colors available: black; blue; red; yellow; or white	18 Spoke (Pair)	\$697.00 <input type="checkbox"/>
Rear Wheels with Hub Locking Brakes- K0108			\$550.00 <input type="checkbox"/>
	Standard pneumatic or solid rubber rear wheels only Cannot be used with Frame Reinforcement		
Rear Wheels with Drum Brakes- K0108			\$1,995.00 <input type="checkbox"/>
	For attendant brake control. Standard pneumatic or solid rubber rear wheels only		


Brake Lever Extension- 2.4"- E0961				
	Simpler accessibility and less strength necessary Standard pneumatic or solid rubber rear wheels only	Right Left	\$75.00 \$75.00	<input type="checkbox"/> <input type="checkbox"/>
Camber- K0108				
	Improved driving performance and agility Cannot be combined with Hub locking Brakes, Drum Brakes, One Arm Drive or E-Motion Wheels	3° 6°	\$274.00 \$274.00	<input type="checkbox"/> <input type="checkbox"/>
Knee Support Options				
Rectangular- K0108		Conical- K0108		
	2" x 4" \$35.00 <input type="checkbox"/> 2.7" x 4" \$35.00 <input type="checkbox"/> 3.1" x 4" \$35.00 <input type="checkbox"/> 2" x 6" \$35.00 <input type="checkbox"/> 3" x 6" \$35.00 <input type="checkbox"/>		2" x 4" x 4" \$44.00 <input type="checkbox"/> 3.5" x 5.5" x 4" \$44.00 <input type="checkbox"/> 5" x 7" x 4" \$44.00 <input type="checkbox"/>	
Footplate Options				
Choose Footplate Angle				
	Angle Adjustable Footplate Set in Position 1- 80°		No charge	<input type="checkbox"/>
	Angle Adjustable Footplate Set in Position 2- 90°		No charge	<input type="checkbox"/>
	Angle Adjustable Footplate Set in Position 3- 100°		No charge	<input type="checkbox"/>
1-Piece Fixed -- No Flip-Up to Rear- K0040			\$140.00	<input type="checkbox"/>
				
1-Piece Rigid -- Flip-Up to Rear- K0040			\$199.00	<input type="checkbox"/>
				
2-Piece -- Dual Flip-Up to Sides- K0040			\$258.00	<input type="checkbox"/>
				
Footrest Accessories				
Seat Depth / Footrest Quick Adjustment- K0108			\$118.00	<input type="checkbox"/>
	Allows for quick adjustment of seat depth and footrest height			
Calf Strap- K0038			\$73.00	<input type="checkbox"/>
Must be ordered				
Footrest Height Adapter- K0108				
	Clients with hip width 13.8"/15.0"/17.7"/19"	Pair	\$63.00	<input type="checkbox"/>
	Clients with hip width 15.7"/16.5"	Pair	\$205.00	<input type="checkbox"/>
Cannot be used with 2-Piece Dual Flip-Up Footrest / Must order 1-Piece No Flip-Up to Rear Footrest				
Removable Legrest / Quick Release- K0108			\$135.00	<input type="checkbox"/>
Allows legrest / footrest assembly to be removed using pins- easier transport				
Foot Guides- K0108				
	Allows individual foot adjustment to accommodate client needs while increasing stability	Small Right Small Left	\$213.00 \$213.00	<input type="checkbox"/> <input type="checkbox"/>
		Large Right Large Left	\$233.00 \$233.00	<input type="checkbox"/> <input type="checkbox"/>

Footrest Wedge- K0108				
	Wedge to help straighten and stabilize feet	2.0" x 4.0" x 3.1"	\$76.00	<input type="checkbox"/>
		3.1" x 5.0" x 3.1"	\$84.00	<input type="checkbox"/>
Non-Slip Footrest Surface- K0108			\$129.00	<input type="checkbox"/>
Accessories				
Foam Seat Cushion with Seat Cushion Cover- E2601				
Additional Seat Cushion Cover- E2619				
	Seat Cushion		\$397.00	<input type="checkbox"/>
	Additional Seat Cushion Cover		\$219.00	<input type="checkbox"/>
Steel Handrims- E0967			\$120.00	<input type="checkbox"/>
				
Maxgrip Handrims- K0108				
	Handrims- rubber covered			
	Original	\$578.00	<input type="checkbox"/>	
	Ergo Para		\$915.00	<input type="checkbox"/>
				
One Arm Drive- E0958				
	Standard wheel only	Right	\$1,595.00	<input type="checkbox"/>
		Left	\$1,595.00	<input type="checkbox"/>
Table Tray / Fixed- E0950			\$378.00	<input type="checkbox"/>
	Requires special backrest. Call for description			
Push Handles- K0108			\$159.00	<input type="checkbox"/>
	Standard Push Handles			
Push Handle Extension- K0108			\$307.00	<input type="checkbox"/>
	Extends Push Handles Adjustable 1" - 6"			
Traditional Anti-Tips - E0971			\$152.00	<input type="checkbox"/>
				
Swing Away Anti-Tips- K0108			\$450.00	<input type="checkbox"/>
				
XL Adapter- Anti-Tip- K0108			\$342.00	<input type="checkbox"/>
Needed for rear wheel position XL when using the Swing Away Anti-Tips Cannot be used with Frame Reinforcement				
Tip Aid for Attendant- K0108			\$57.00	<input type="checkbox"/>
	Tip aid is for attending persons to help to lift up the front wheels climbing a curb			

	<p>Spoke Guard-Transparent- K0065</p> <p>Plexiglass shields to protect hands and fingers</p>	<p>\$186.00 <input type="checkbox"/></p>
-----------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------	-------------------------------------------------

Outdoor / Travel Accessories

	<p>Vehicle Tie-Down Kit- K0108</p> <p>The Tie-Down Kit must be used for client transportation in the car. Chair has to be tied down using the Q'Straint 5001-T2 or any restraint system certified by Q'Straint 5001-T2 . The occupant has to be safely secured using the Q'Straint Vehicle Anch 3 point or any restraint system certified by ISO 10542/2 standards. Frame Reinforcement and Swing Away Anti-Tips are REQUIRED</p>	<p>\$988.00 <input type="checkbox"/></p>
-----------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------

	<p>Frame Reinforcement- K0108</p> <p>Frame Reinforcement is needed to attach the Vehicle Tie Down Kit or E-Motion Power Assist system. Must order the Swing Away Anti-Tips Cannot be used with Hub Locking Brakes or Traditional Anti-Tips</p>	<p>\$290.00 <input type="checkbox"/></p>
-----------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------

Chair Order Policy

Contact LEVO USA at 1-888-LEVOUSA with any questions regarding how to accurately complete this form.

LEVO USA is not responsible for configuration or size discrepancies resulting from inaccurate or incomplete information provided on the order form.

Order confirmations will be sent to you. Please review and sign your agreement/confirmation authorizing submission of the order as stated. Return the signed copy to LEVO USA at fax number 763-582-0442 or via email at: amyj@danetechnologies.com.

All chairs are built to specifications provided, therefore, all sales are final. No returns or exchanges. Your signature acknowledges your understanding and acceptance of this no return policy. Customer is responsible for full payment of chair shipped based on specifications provided on confirmed order form.

The HCPC codes identified herein are our best interpretation of the match between chair, accessories and the appropriate code. This coding interpretation has not been verified by DME PDAC nor should our interpretation replace the customer's own responsibility for submitting accurate information to third party sources.

Please note: All prices and specifications provided on form are subject to change without prior notice. All prices for options or accessories on order form are only valid if ordered with the chair or prior to completed production of unit. Any chair modifications, options or accessories required for future use or maintenance of chair will carry different pricing and warranty terms.

Warranty Policy

Your LEVO product is guaranteed from the date of purchase for:

Two years covering all material and manufacturing defects of mechanical parts.

LEVO USA will not repair or replace free of charge any part or parts found to be defective due to abuse, misuse or lack of maintenance.

