



Client Name: _____

Date: _____

Height: _____ ft / in

Diagnosis

- Paraplegia Quadriplegia Hemiplegia
 Multiple Sclerosis Muscl. Dystrophy Cerebral Palsy
 Other:

Weight: _____ lbs (265lb max)

Reasons for purchasing a LEVO

Practical / Functional

- Only for standing Use as primary wheelchair
 Use at home Use as secondary wheelchair
 Use at work / school Other:

Medical / Therapeutical

- Decubitus Contractures Spasms
 Digestion Bladder problems Osteoporosis
 Circulation Pain (postural pain) Other:

How did you first hear of LEVO?

- Dealer Therapist Internet
 Advertisement Doctor
 Other: Friends

Financing

- Auto / Disability Medicare Work Comp
 Medicaid / State Private Insurance Donations
 Personal Voc. Rehab

Dealer Information

Contact _____

State _____

Dealers Name _____

Fax# _____

Address _____

Ph# _____

City _____

E-mail _____

Dealer Code _____

P.O.# _____

LCEV/LCEV KID & JR Base

Includes: Power stand and anti puncture solid rubber front casters; manual drive base; safety side guards; extra strong knee support; basic frame color; armrests round with built-in control switch, rear air filled tires with aluminum hand rims and standard locking brakes. Crash tested per ISO 7176/19 ANSI/RESNA WC/Vol. 1-Section 19

Price **\$12,875.00**

HCPC


Batteries- Required

\$130.00


Basic frame colors

- Platinum **No charge** Red **No charge** Black **No charge**
 Yellow **No charge** Blue **No charge** Chrome **\$950.00**

Please answer the following questions accurately and completely. Identifying the exact measurements within the diagrams allows for correct configuration and final adjustments for the client. Please do not provide measurements from client's existing chair or cushion!



- B. _____ Upper leg length
 C. _____ Lower leg length
 D. _____ Foot length
 F. _____ Seat to top of shoulder
 H. _____ Seat to top of head
 I. _____ Arm height
 K. _____ Chest depth



- A. _____ Hip width
 E. _____ Seat to scapula
 G. _____ Shoulder width


Non-Standard Seat Frame Measurements

- If "A" above measures 19"
 If "B" above measures 20"-22"
 If "B" above measures 22"-24"
 If "C" above measures 19"-22" Requires 26" rear wheels w/6" front casters
 If "A" above measures 14" on the LCEV KID/JR chair only

- \$1,400.00**
\$218.00 **E2203**
\$218.00 **E2204**
\$218.00
\$1,400.00

Backrest Options

Backrest height



Standard-

Measure the distance between the seating surface and the lower shoulder blade bone

R = _____ inches

Omit Backrest option

Tension Adjust Sling back Backrest -----	No charge	<input type="checkbox"/>	
V-Trak Standard Backrest Low, 16" high -----	\$175.00	<input type="checkbox"/>	
V-Trak Standard, Backrest, 24" high -----	\$1,063.00	<input type="checkbox"/>	E2620
V-Trak Standard, Backrest, 24" high -----	\$1,083.00	<input type="checkbox"/>	E2620
V-Trak High Backrest, 31.5" high -----	\$1,158.00	<input type="checkbox"/>	E2620
V-Trak Axxis Reduced Height Backrest, 12" high -----	\$1,279.00	<input type="checkbox"/>	E2620
V-Trak Axxis Low Backrest, 16" high -----	\$1,304.00	<input type="checkbox"/>	E2620
V-Trak Axxis Standard, Backrest 24" high -----	\$1,550.00	<input type="checkbox"/>	E2620
Additional V-Trak arms -----	\$772.00	<input type="checkbox"/>	
V-Trak Backrest extension 8" (No Upholstery) -----	\$526.00	<input type="checkbox"/>	

V-Trak Headrest (incl. Track extension)

V-Trak Headrest winged (V-Trak backrest-Required)	\$788.00	<input type="checkbox"/>	E0955
V-Trak Headrest simple (V-Trak backrest-Required)	\$730.00	<input type="checkbox"/>	E0955

V-Trak Thoracic support & pad

V-Trak Trunk support & pad & cover L/XL right & left (V-Trak backrest-Required)	\$1,410.00	<input type="checkbox"/>	E0956
V-Trak Trunk support & pad & cover S/M right & left (V-Trak backrest-Required)	\$1,410.00	<input type="checkbox"/>	E0956

Standard Chest belt-Required

Body Point padded chest belt	\$75.00	<input type="checkbox"/>	E0960
------------------------------	---------	--------------------------	-------

Standing Operation

	Right	<input type="checkbox"/>	Left	<input type="checkbox"/>	
Armrest round with hand control					No charge <input type="checkbox"/>
Standard - armrest round with control switch					No charge <input type="checkbox"/>
Armrest 6 cm/2.4" wide with hand control					\$360.00 <input type="checkbox"/>
Armrest 6 cm/2.4" wide with built-in control switch					\$360.00 <input type="checkbox"/>
Armrest 5 cm/2" longer					\$321.00 <input type="checkbox"/>

Swing away chest support arms

Rotate and swing away to the rear. Increased security for the user in standing position	10.5"	\$697.00	<input type="checkbox"/>
Distance between backrest tube and the arm mid-point	12.5"	\$697.00	<input type="checkbox"/>
	14.5"	\$697.00	<input type="checkbox"/>

Dual Switch for table control

Stand-up operation with push buttons for table. Table tray not included (see accessories).	\$1,160.00
--	------------

Control switch extension

Vertical for round armrest	\$102.00	<input type="checkbox"/>
Vertical for wide armrest	\$142.82	<input type="checkbox"/>
Horizontal	\$30.00	<input type="checkbox"/>


Rear Wheel Options & Accessories

Rear wheel size

Rear wheel 20", castor 5" - Kid/JR only	f 13.0" r 18.9"	No charge	<input type="checkbox"/>
Rear wheel 22", castor 5"	f 19.0" r 17.7"	No charge	<input type="checkbox"/>
Rear wheel 24", castor 5"	f 20.0" r 18.5"	No charge	<input type="checkbox"/>
Rear wheel 24", castor 5"	f 21.2" r 19.3"	No charge	<input type="checkbox"/>
Rear wheel 26", castor 6"	f 22.5" r 20.5"	\$295.00	<input type="checkbox"/>

Rear wheel position

Position S Position M Position L Position XL



The shorter the wheel base, the more dynamic the handling of the wheelchair, but also the tilt risk is increased. For example, with wheelchair size "L" you could choose the rear wheel position "M". This combination creates quicker response and the need for anti-tippers.

Rear wheels/solid rubber tires

Aluminum hand ring included	\$145.00	<input type="checkbox"/>	E2220
-----------------------------	----------	--------------------------	-------

Rear wheels with hub locking brakes

Standard pneumatic or solid rubber rear wheels only. Can not be used with Frame Reinforcement	\$550.00	<input type="checkbox"/>
---	----------	--------------------------

Rear wheels with drum breaks

For attendant brake control. Standard pneumatic or solid rubber rear wheels only	\$1,995.00	<input type="checkbox"/>
--	------------	--------------------------

Brake lever extension (6 cm/2.4")

Simpler accessibility and less strength necessary.	Right	\$240.00	<input type="checkbox"/>	E0961
Standard pneumatic or solid rubber rear wheels only.	Left	\$240.00	<input type="checkbox"/>	E0961

Rear Wheel Options & Accessories

Rear wheels "Spinergy"

Light weight, improved driving performance, easy to clean
Includes mounted aluminum hand rim.
Can not be used with Hub Locking brakes or E-Motion wheels
Rim color is black. Spoke colors available: black, blue, red, yellow & white

12 SPOKE (Pair) **\$925.00**
18 SPOKE (Pair) **\$697.00**

Camber

Improved driving performance and agility
Cannot be combined with Hub locking brakes, Drum brakes, One arm drive or E-Motion wheels

3° **\$274.00**
6° **\$274.00**

Leg Rest Option

Knee wedge

To be used with standard knee support provided in base price. Supports leg angle

Rectangular

2.0" x 4" **\$25.00**
2.7" x 4" **\$25.00**
3.1" x 4" **\$25.00**
2.0" x 6" **\$25.00**

Conical- More support, especially with spasticity

2.0" - 4.0" x 4" **\$44.00**
3.5" - 5.5" x 4" **\$44.00**
5.0" - 7.0" x 4" **\$44.00**

Footrest Options

All footrest options are adjustable- Choose footrest angle

80° 90° 100° No charge **K0040**

One-piece footrest fix, no flip-up feature **\$195.00**

Two-piece footrest with flip-up feature **\$195.00**

Ground level footrest/Large size **\$239.00** **K0041**

Seat depth / Footrest quick adjustment **\$117.00**

Calf Strap-Required **\$32.00** **K0038**

Footrest height adapter

Clients with hip width 13.8"/15.0"/17.7/19"
Clients with hip width 15.7"/16.5"- Cannot be used with 2-piece flip-up footrest/Must order 1-piece, No flip footrests

(Pair) **\$62.00**
(Pair) **\$202.00**

Removable Leg rest/Quick release **\$133.00**

Non-slip footrest surface **\$129.00**

Foot guide

Can be adjusted individually. Secures and aligns feet

Small right **\$213.00** Large right **\$233.00**
Small left **\$213.00** Large Left **\$233.00**

Knee support height adapter

Is suitable for persons with long lower legs. Knee holder gets closer to the joint
Better pressure distribution

Pair **\$298.00**
Each R or L **\$149.00**

Accessories

Foam seat cushion with cover

Omit seat cushion option ----- **No charge**
2" Foam seat cushion with cover----- **\$277.00**
Additional seat cushion cover ----- **\$128.00** **E2619**

Backrest Extension - 6"

Only used with Standard Tension Adjust Sling Back backrest - includes upholstery

\$331.00

Back rest extension with push handle and upholstery

Length 5.9" ----- **\$739.00**
Length 9.8" ----- **\$755.00**
Length 13.8" ----- **\$778.00**

Accessories

Stroller handle/Push Handle Bar

Can be used with standard or V-Track backrest

\$297.00

Stroller handle Height adjustable for Attendant Person - Kid/Jr only

Handle length 12" - 14"

\$350.00

Stroller handle Height adjustable with up/down switch for Attendant- Kid/Jr only

Handle length 12" - 14"

\$618.00

Stroller handle - Folding down; Width adjustable - Kid/Jr only

\$826.00

Table tray

Fixed tube in armrest

\$560.00 E0950

Swing away table tray

Support installed at the armrest. The table is still fastened on the chair, even when not in use
Only for Combi arm rest 2.4" / 6.1 cm width

\$1,102.00 E0950

Steel hand ring

Increased grip, less temperature sensitive, cleaner

\$119.00 E0967

Maxgrip

Hand ring-rubber covered

Original (Pair)

\$570.00

Ergo Para (Pair)

\$905.00

One arm drive

Standard wheel only

Right
Left

\$1,916.00 E0958
\$1,916.00 E0958

Safety guards

Right & Left (Pair)

No charge

Push handles

Standard push handles

\$239.00

Push handle extension adjustable

Extends push handles

\$595.00

Traditional Anti-tip rolls

Right & Left (Pair)

\$266.00 E0971

Swing away Anti-tip "HC"

Must be used with Frame Reinforcement option

Right & Left (Pair)

\$450.00 E0971

XL adapter- Anti-Tip

Needed for rear wheel position XL when using the Swing away Anti-tip "HC"

Right & Left (Pair)

\$342.00 E0971

Tip aid for Attendant

Tip aid is for attending persons to help to lift up the front wheels climbing a curb

\$57.00 E0971

Spoke guard-Transparent

Plexi glass shields to protect hands and fingers

Right & Left (Pair)

\$186.00 K0065

Monitor

Provides an objective outcome measurement for users and therapists. Delivered with PC-Software

\$690.00

Vehicle Mount & Driving Accessories

Vehicle Tie Down Kit

The tie-down kit must be used for client transportation in the car. Chair has to be tied down using the Q'Straint 5001-T2 or any restraint system certified by ISO 10542/2 standards.

The occupant has to be safely secured using the Q'Straint Vehicle Anch 3 point or any restraint system certified by ISO 10542/2 standards

\$972.00

Frame Reinforcement

Frame reinforcement is needed to attach the E-Motion Power Assist system. Must order the swing away "HC" anti-tippers when using the E-Motion Power Assist system

\$270.00

Chair Order Policy

Contact LEVO USA at 1-888-LEVOUSA with any questions regarding completion of this form.

LEVO USA is not responsible for configuration or size discrepancies resulting from inaccurate or incomplete information provided on the order form.

Order confirmations will be sent to you. Please review and sign your agreement/confirmation authorizing submission of the order as stated. Return the signed copy to LEVO USA at fax number 763-582-0442 or via email at: amyj@danetechnologies.com

All sales are final. No returns or exchanges. Customer is responsible for full payment of chair shipped based on specifications provided on confirmed order form.

The HCPC codes identified herein are our best interpretation of the match between chair, accessories and the appropriate code. It is not a legal determination nor should our interpretation replace customer's own responsibility for submitting accurate information to third party sources.

Please note: All prices and specifications provided on form are subject to change without prior notice. All prices for options or accessories on order form are only valid if ordered with the chair or prior to completed production of unit. Any chair modifications, options or accessories required for future use or maintenance of chair will carry different pricing and warranty terms.

Warranty Policy

Your LEVO product is guaranteed from the date of purchase for:

Two years covering all material and manufacturing defects of mechanical parts.

LEVO USA will not repair or replace free of charge any part or parts found to be defective due to abuse, misuse or lack of maintenance.

LEVO Stands for a better quality of life



LEVO ^{dt} 
THE EXPERTS IN STANDING

