

Standing Operation

Manual Wheelchair Accessory, **MANUAL** Standing System

E2230 \$8,635



Wheelchair Accessory, **POWER** Standing System, Any Type

E2301 \$10,795



Client weight

(This section applies to MANUAL Standing option only)

Client weight:

Gas strut selection chart

| Depth: ▶ Weight: ▼ | 11.5" 29.0 cm | 12" 31.2 cm | 13" 33.4 cm | 14" 35.6 cm | 15" 37.8 cm | 16" 40.0 cm | 17" 43.0 cm | 18" 45.2 cm | 18.5" 47.4 cm | 19.5" 49.6 cm | 20" 51.8 cm |
|-----------------------|------------------|----------------|----------------|----------------|----------------|----------------|----------------|----------------|------------------|------------------|----------------|
| 88 lbs. to 132 lbs. | 400n | 400n | 400n | 400n | 400n | 400n | 400n | 400n | 400n | 400n | |
| 132 lbs. to 176 lbs. | | | 400n | 500n | 500n | 500n | 500n | 500n | 500n | 500n | |
| 176 lbs. to 220 lbs. | | | | 500n | 600n | 600n | 600n | 600n | 600n | 600n | 600n |
| 220 lbs. to 264 lbs. | | | | | | | 600n | 700n | 700n | 700n | 700n |

Frame color



Black
Blue
Champagne
Red
Yellow

| | | |
|----------------------|----------------------|---|
| <input type="text"/> | <input type="text"/> | • |
| <input type="text"/> | <input type="text"/> | • |
| <input type="text"/> | <input type="text"/> | • |
| <input type="text"/> | <input type="text"/> | • |



Blue metallic
Green metallic
Red metallic
Silver metallic

| | | |
|--------|----------------------|--------|
| Option | <input type="text"/> | \$ 595 |
| Option | <input type="text"/> | \$ 595 |
| Option | <input type="text"/> | \$ 595 |
| Option | <input type="text"/> | \$ 595 |

Orange reflex
Lemon reflex
Silver reflex

| | | |
|--------|----------------------|--------|
| Option | <input type="text"/> | \$ 695 |
| Option | <input type="text"/> | \$ 695 |
| Option | <input type="text"/> | \$ 695 |

Aluminium polished
Traffic white
Other RAL color

RAL: 9016

| | | |
|--------|----------------------|--------|
| Option | <input type="text"/> | \$ 995 |
| Option | <input type="text"/> | \$ 895 |
| Option | <input type="text"/> | \$ 895 |

Standing Release

(This section applies to MANUAL Standing option only)

| | | |
|------------------|----------------------|---|
| Bilateral | <input type="text"/> | • |
| Right | <input type="text"/> | • |
| Left | <input type="text"/> | • |

Armrest

(This section applies to Power Standing System only)


Armrest - Angle adjustable

Arm height

Note: Please consider cushion thickness when choosing the arm height

Seat to elbow:


5.5" min Set to: 8.75" max



- 1 Armrest - long (15" from backpost), round with foam pad, height adj.
- 2 Armrest - short (13" from backpost), round with foam pad, height adj.
- 3 Armrest - long 2.4" wide polyurethane pad, height adj.
- 4 Armrest - short 2.4" wide polyurethane pad, height adj.
- 5 Armrest - Swing-away chest support, with arms, height adj.

NONE

| | | | |
|--------|-------|----------------------|-------|
| Option | E0973 | <input type="text"/> | \$107 |
| Std. | E0973 | <input type="text"/> | • |
| Option | E0973 | <input type="text"/> | \$366 |
| Option | E0973 | <input type="text"/> | \$251 |
| Option | E0973 | <input type="text"/> | \$753 |



Flip up and rotate out, with round tube and foam pad only

Distance to beginning of bend, see picture at left

Choose "D" dimension:
10.5"
12.5"
14.5"



Standing control

- Switch in arm - left
- Switch in arm - right
- Pendant hand control

(This section applies to Power Standing System only)

| | | |
|--------|----------------------|---|
| Option | <input type="text"/> | • |
| Option | <input type="text"/> | • |
| Std. | <input type="text"/> | • |



Batteries

Standard Battery Box

K0108 \$172

7

Lithium (Wheelchair is approximately 6 lbs. lighter - Not for use with client weight 185 lbs or higher in combination with 19" or 19.75" wide chairs or in combination with seat depths 19.5"; 20" or 21")

K0108 \$730

Backrest


Backrest - height adjustable canes with tension adjustable sling and cushion

| Upholstery height: | | |
|--------------------|--|---|
| 12.25" | | • |
| 13.00" | | • |
| 13.75" | | • |
| 14.5" | | • |
| 15.5" | | • |
| 16" | | • |

Back - V-Trak "Dynamic Stability" with internal adjustable straps & laterals (flexes to absorb energy, .25" padded/breathable up.)
(must order mounting hardware)

V-TRAK Axis with parallel sides

Parallel Sides with 2-3 inches of contour depth




Internal adjustable spine profiling straps



Internal adjustable lateral support

| Width | For Seat Width | | For Seat Widths | | For Seat Widths | |
|-----------------|----------------|-------|--------------------------|-------|----------------------------|-------|
| | 13.75" (M) | | 15" / 15.75" / 16.5" (L) | | 17.75" / 19" / 19.75" (XL) | |
| 1 Height 12" | \$695 | E2620 | \$695 | E2620 | \$750 | E2620 |
| 2 Height 15.75" | \$695 | E2620 | \$695 | E2620 | \$750 | E2620 |
| 3 Height 21.5" | \$750 | E2620 | \$775 | E2620 | \$795 | E2620 |

Note: The 21.5" high backrest can't be used on seat width of 13.75"



V-TRAK Axis Deep Contour (Lateral crest is always 22cm from bottom)

Deep contour with 5-7 Inches of contour depth



Internal adjustable spine profiling straps


Internal adjustable lateral support - Lateral adjustment for front of backrest

| Width | For Seat Width | | For Seat Widths | | For Seat Width | | For Seat Width | |
|-----------------|----------------|-------|--------------------------|-------|----------------|-------|--------------------|-------|
| | 13.75" (M) | | 15" / 15.75" / 16.5" (L) | | 17.75" (XL) | | 19" / 19.75" (XXL) | |
| 4 Height 12" | \$695 | E2620 | \$695 | E2620 | \$750 | E2620 | \$850 | E2620 |
| 5 Height 15.75" | \$775 | E2620 | \$795 | E2620 | \$825 | E2620 | \$850 | E2620 |
| 5 Height 21.5" | \$775 | E2620 | \$795 | E2620 | \$825 | E2620 | \$850 | E2620 |

Note: The 21.5" high backrest can't be used on seat width of 13.75"



Mounting Hardware - Adjustable



| | | | |
|---|------|-------|--------|
| 1 Lite - tube style | NONE | E1028 | \$ 292 |
| 2 Original - flat style (recommended for chair back widths over 19") | NONE | E1028 | \$ 292 |

Backrest options




| | | | |
|--|------|-------|--------|
| 3 Standard Backrest - 6" tube extension for sling back (for raising back height 17" to 22") | NONE | K0108 | \$ 369 |
| V-Trak backrest - 8" extension | | K0108 | \$ 543 |

Backrest angle

| | | |
|-----|--|---|
| 85° | | • |
| 90° | | • |
| 95° | | • |

Seat Cushion - memory foam with cover



| Seat width | 13.75" | 15" | 15.75" | 16.5" | 17.75" | 19" |
|---|--------|-------|--------|-------|--------|-------|
| | 35 cm | 36 cm | 37 cm | 38 cm | 39 cm | 40 cm |
| Seat depth | | | | | | |
| Matches seat depth unless otherwise specified | | | | | | |

E2601 \$ 429

Additional Seat Cushion Cover

| | | | |
|-------------|------|-------|--------|
| Standard | NONE | E2619 | \$ 237 |
| Incontinent | NONE | E2619 | \$ 237 |

Footplate angle

Set to:

| | | |
|------|-------|--|
| 80° | K0040 | |
| 90° | K0040 | |
| 100° | K0040 | |

Accessories



| | | | | | |
|----|---|------|--|-------|----------|
| 1 | Tray table <i>(applies to POWER Standing System only)</i> | | | | |
| | Tray table - hardware | | | E0950 | \$ 305 |
| | | | | E1028 | \$ 316 |
| | Push handles | Pair | | K0108 | \$ 189 |
| 2 | Push handles - 2.75" extension | Pair | | K0108 | \$ 365 |
| 3 | Push handles - detachable, height adjustable 1" - 6" | Pair | | K0108 | \$ 378 |
| 4 | Anti Tips - traditional- Pair | Pair | | E0971 | \$ 164 |
| 5 | Tip Aid - attendant- Pair | Pair | | K0108 | \$ 92 |
| 6 | Anti Tips - swing away with attendant tip aid | Pair | | K0108 | \$ 486 |
| | XL Adapter* <i>"Necessary with swing away anti-tips for seat depth positions- of 19.5", 20" & 21" when not using E-Fix or E-Motion wheels- Use of E-Fix or E-Motion wheels in these seat depths require Anti-tips provided by Alber</i> | Pair | | K0108 | \$ 164 |
| 7 | Push to Lock- Wheel lock extensions - 2.5" | Pair | | K0961 | \$ 162 |
| | | | | | |
| 8 | Seat depth and footrest quick adjustment kit <i>Includes 8 clamps with post and 4 threaded knobs</i> | | | K0108 | \$ 139 |
| 10 | Frame reinforcement kit <i>Necessary when e-fix or similar system is mounted Anti-tips - swing away required</i> | | | K0108 | \$ 313 |
| 11 | Spoke Guards- Transparent | Pair | | K0065 | \$ 201 |
| | Vehicle tie-down kit | | | K0108 | \$ 1,067 |

*Fixed attachment points for chair transport Includes frame reinforcement kit
Anti-tips - swing away required- Note: When using E-Fix or E-Motion wheels for seat depth positions of 19.5", 20" & 21" Anti-tips provided by Alber are required*

Mandatory if client is transported in vehicle. Chair required to be secured with the Q'Straint 5001-T2 or any restraint system certified by Q'Straint 5001-T2. Occupant must be secured using the Q'Straint Vehicle 3 point or any restraint system certified by ISO 10542/2 standards.

Remarks