



Dealer: _____ Contact: _____ Date: _____
 Address/City/St/Zip: _____
 Phone: _____ FAX: _____ Email: _____
 Clinician: _____ Phone: _____ Email: _____
 Client: _____ Disability: _____ Misc.: _____

LEVO USA
 10400 Fountains Drive
 Maple Grove, MN 55369
 Toll Free: 888-538-6872
 FAX: 763-582-0442
 Orders to:
 sales@levousa.com


Wheelchair Base K0005 \$3,235

- Manually Propelled Ultralight Wheelchair
- Maximum Client Weight: 265 lbs.
- Tension Adjustable Seat & Back Upholstery
- Fold Down, Angle Adjustable Backrest, Height Adjustable 12.5"-16"
- Rear Wheel - 24", high pressure pneumatic, hub black, 18 spoke silver, w/aluminum hand rim
- Front caster - solid rubber, composite wheel
- Standard colors: platinum, red, yellow, blue or black
- 1 Piece Footplate - Fixed
- 0 Degree Rear Wheel Camber
- Standard Push to Lock Brakes
- Side Guards



Anatomical Dimensions

- Weight: _____
- Height: _____
- Hip Width: ①
- Upper Leg Length: ②
- Lower Leg Length: ③
- Seat to Elbow: ④
- Chest Width: ⑤
- Chest Depth: ⑥
- Floor to top of knee: ⑦
- Seat to Top of Shoulder: ⑧



① Hip width > Seat width

Hip width: _____

①
Seat width: _____

35 cm	13.75"				
38 cm	15"				
40 cm	15.75"				
42 cm	16.5"				
45 cm	17.75"				
48 cm	19"				\$ 1,035
50 cm	19.75"				\$ 1,035

(Upcharge for 19" seat width applies to **Manual Standing System**)
 (19.75" width available only on **Power Standing System chair**)

② Upper leg length/Seat depth (upholstery)

With proper seat depth, the ball of the knee should be vertically aligned with the Knee Block Mounting Post

Upper leg length: _____

The seat depth of the upholstery is approximately 1.5" to 2" from the back of the knee

Seat depth (upholstery): _____

Holes showing at rear of seat rail

29.0 cm	11.5"	0			
31.2 cm	12"	1			
33.4 cm	13"	2			
35.6 cm	14"	3			
37.8 cm	15"	4			
40.0 cm	16"	5			
43.0 cm	17"	0			
45.2 cm	18"	1			
47.4 cm	18.5"	2			
49.6 cm	19.5"	3			E2203
51.8 cm	20"	4			E2203
54.0 cm	21"	5			

There are two ranges seat of depth adjustment for the LAE (S-M & L-XL)
 Depth adjustment is made utilizing the telescoping upper and lower seat rails

③ Lower leg length/Seat to footplate

The *Lower leg length* determines the *Seat to footplate* setting, possible wheel/caster combinations and front seat heights.

Lower leg length: _____

Please take the cushion into consideration when choosing the *Seat to footplate* setting.

Seat Height Front:	19"	20"	21"	22.5"
Rear Wheel	22"	24"	24"	26"
Seat to footplate:	5" caster	5" caster	5" caster	6" caster
14"				
15"				
15.75"				
16.5"				
17.25"				
18"				
19"				
19.75"				
20.5"				
21.25"				

Rear seat height is approximately 1.5" lower than front

Mandatory with 19" seat height
 Tire - 22" x 1" rear wheel- hub black, spoke silver with 5" front caster

Pair \$ -

Mandatory with 22.5" seat height
 Tire - 26" x 1" rear wheel- hub black, spoke silver with 6" front caster

NONE

Pair \$ 351

Standing Operation

Manual Wheelchair Accessory, **MANUAL** Standing System

E2230 \$8,635



Wheelchair Accessory, **POWER** Standing System, Any Type

E2301 \$10,795



Client weight

(This section applies to MANUAL Standing option only)

Client weight:

Gas strut selection chart

Depth: ▶ Weight: ▼	11.5" 29.0 cm	12" 31.2 cm	13" 33.4 cm	14" 35.6 cm	15" 37.8 cm	16" 40.0 cm	17" 43.0 cm	18" 45.2 cm	18.5" 47.4 cm	19.5" 49.6 cm	20" 51.8 cm
88 lbs. to 132 lbs.	400n	400n	400n	400n	400n	400n	400n	400n	400n	400n	
132 lbs. to 176 lbs.			400n	500n	500n	500n	500n	500n	500n	500n	
176 lbs. to 220 lbs.				500n	600n	600n	600n	600n	600n	600n	600n
220 lbs. to 264 lbs.							600n	700n	700n	700n	700n

Frame color

All color choices marked with • are no charge



Black
Blue
Wine Red
Red
Yellow

<input type="checkbox"/>		•
<input type="checkbox"/>		•
<input type="checkbox"/>		•
<input type="checkbox"/>		•
<input type="checkbox"/>		•



Blue metallic
Basalt Grey
Red metallic
Silver metallic

Orange reflex
Lemon reflex
Silver reflex

White Aluminum
Traffic white
Other RAL color

<input type="checkbox"/>	Option		\$	495
<input type="checkbox"/>	Option		•	
<input type="checkbox"/>	Option		\$	495
<input type="checkbox"/>	Option		\$	495
<input type="checkbox"/>	Option		\$	595
<input type="checkbox"/>	Option		\$	595
<input type="checkbox"/>	Option		\$	595
<input type="checkbox"/>	Option		\$	•
<input type="checkbox"/>	Option		•	
<input type="checkbox"/>	Option		\$	695

Standing Release

(This section applies to MANUAL Standing option only)

Bilateral- Recommended

<input type="checkbox"/>		•
<input type="checkbox"/>	Right	•
<input type="checkbox"/>	Left	•

Armrest

(This section applies to Power Standing System only)


Armrest - Angle adjustable

Arm height

Note: Please consider cushion thickness when choosing the arm height

Seat to elbow:

5.5" min Set to: 8.75" max



- 1 Armrest - long (15" from backpost), round with foam pad, height adj.
- 2 Armrest - short (13" from backpost), round with foam pad, height adj.
- 3 Armrest - long 2.4" wide polyurethane pad, height adj.
- 4 Armrest - short 2.4" wide polyurethane pad, height adj.
- 5 Armrest - Swing-away chest support, with arms, height adj.


NONE

<input type="checkbox"/>	Option	E0973		\$107
<input type="checkbox"/>	Std.	E0973		•
<input type="checkbox"/>	Option	E0973		\$366
<input type="checkbox"/>	Option	E0973		\$251
<input type="checkbox"/>	Option	E0973		\$753

Flip up and rotate out, with round tube and foam pad only

Distance to beginning of bend, see picture at left

Choose "D" dimension:
10.5"
12.5"
14.5"



Standing control

- Switch in arm - left
- Switch in arm - right
- Pendant hand control

(This section applies to Power Standing System only)

<input type="checkbox"/>	Option		•
<input type="checkbox"/>	Option		•
<input type="checkbox"/>	Std.		•

Batteries

Standard Battery Box

K0108 \$172

Lithium (Wheelchair is approximately 6 lbs. lighter)- Not for use with client weight 185 lbs or higher or for 48cm or 50cm wide chairs

K0108 \$730



Backrest

Backrest - height adjustable canes with tension adjustable sling and cushion

Upholstery height:		
12.25"		•
13.00"		•
13.75"		•
14.5"		•
15.5"		•
16"		•

Back - V-Trak "Dynamic Stability" with internal adjustable straps & laterals (flexes to absorb energy, .25" padded/breathable up.)
(must order mounting hardware)

V-TRAK Axis with parallel sides

Parallel Sides with 2-3 inches of contour depth




Internal adjustable spine profiling straps



Internal adjustable lateral support

Width	For Seat Width		For Seat Widths		For Seat Widths	
	13.75" (M)		15" / 15.75" / 16.5" (L)		17.75" / 19" / 19.75" (XL)	
1 Height 12"	\$695	E2620	\$695	E2620	\$750	E2620
2 Height 15.75"	\$695	E2620	\$695	E2620	\$750	E2620
3 Height 21.5"	\$750	E2620	\$775	E2620	\$795	E2620

Note: The 21.5" high backrest can't be used on seat width of 13.75"



V-TRAK Axis Deep Contour (Lateral crest is always 22cm from bottom)

Deep contour with 5-7 Inches of contour depth



Internal adjustable spine profiling straps


Internal adjustable lateral support - Lateral adjustment for front of backrest

Width	For Seat Width		For Seat Widths		For Seat Width		For Seat Width	
	13.75" (M)		15" / 15.75" / 16.5" (L)		17.75" (XL)		19" / 19.75" (XXL)	
4 Height 12"	\$695	E2620	\$695	E2620	\$750	E2620	\$850	E2620
5 Height 15.75"	\$775	E2620	\$795	E2620	\$825	E2620	\$850	E2620
5 Height 21.5"	\$775	E2620	\$795	E2620	\$825	E2620	\$850	E2620

Note: The 21.5" high backrest can't be used on seat width of 13.75"



Mounting Hardware - Adjustable



1 Lite - tube style	NONE	E1028	\$ 292
2 Original - flat style (recommended for chair back widths over 19")	NONE	E1028	\$ 292

Backrest options




3 Standard Backrest - 6" tube extension for sling back (for raising back height 17" to 22")	NONE	K0108	\$ 369
V-Trak backrest - 8" extension		K0108	\$ 543

Backrest angle

85°		•
90°		•
95°		•

Seat Cushion - memory foam with cover



Seat width	13.75"	15"	15.75"	16.5"	17.75"	19"
	35 cm	36 cm	37 cm	38 cm	39 cm	40 cm
Seat depth						
Matches seat depth unless otherwise specified						

E2601 \$ 429

Additional Seat Cushion Cover

Standard	NONE	E2619	\$ 237
Incontinent		E2619	\$ 237

Supports




1	Knee support - Padded and reinforced								K0108		\$	186
2	Rectangular knee wedge										\$	38
											\$	38
											\$	38
											\$	38
											\$	38
											\$	38
											\$	38
											\$	38
											\$	38
											\$	38
											\$	38
											\$	38
4	Knee support - Sport								K0108		\$	319
5	Knee support - Polyurethane pad, rigid tube, swingaway										\$	589
	Gelovations knee block pad mount with cover for-swing away										\$	589
	Calf strap (shown in pict #1 below)										\$	325
	Chest support belt										\$	79
6	Additional Chest support - BodyPoint Universal Neoprene- 4.5" wide - SP110L									X	\$	130
										X	\$	144
											\$	130
											\$	144
7	Shoe Holder	3.75" x 7.25"	Small	Right					K0108		\$	230
		3.75" x 7.25"	Small	Left					K0108		\$	230
		4.5" x 8"	Large	Right					K0108		\$	252
		4.5" x 8"	Large	Left					K0108		\$	252
	Hip Belt - Airplane pull								E0978		\$	129

Footrest



1	One piece footplate, fixed, angle adjustable								Std.		•	
2	Two-piece footplate, flip-up, angle adjustable								K0108		\$	308
3	One piece footplate, flip-up (back), angle adjustable								K0040		\$	237
2	Footrest height adapter <i>(footplate will be 1" narrower than normal)</i>										\$	75
											\$	243
4	Easy transfer footplate <i>(extended "L" brackets bring footplate closer to the floor)</i>										\$	419
											\$	419
	Detachable legrest/downposts								K0108		\$	161
	Non-slip footplate surface								K0108		\$	139

Wheel options



			Tire				Handrim			Spoke Color						
	Wheel Type		1" Gray High Pressure Pneumatic 1" Black Pneumatic Tires \$145	1" Black Spiny Solid Poly \$172	1" BLK Marathon Solid Poly \$172	Steel \$132	Steel - embossed \$578	Aluminum Clear Anodized	Maxgrepp - Rubber Coated \$578	Maxgrepp - Ergo Para CP \$995	Natural Fit \$578	Natural Fit /Thumbguard \$995	Natural Fit SLT \$578	White	Red	Black
Rear Wheel	22"	Standard 18 Spoke - 22", hub black, spoke silver														
		Spinyer 12 Spoke - 22", hub silver \$925														
		Spinyer 18 Spoke - 22", hub silver \$697														
24"	Standard 18 Spoke - 24", hub black, spoke silver															
	Spinyer 12 Spoke - 24", hub silver \$925															
	Spinyer 18 Spoke - 24", hub silver \$697															
26"	Standard 18 Spoke - 26", hub black, spoke silver															
	Spinyer 12 Spoke - 26", hub silver \$925															
	Spinyer 18 Spoke - 26", hub silver \$697															
4	Camber	0°	3°	6°	Pair								K0108		\$	296
5	Frog Legs-Black Carbon Fiber Forks & Casters												K0108		\$	525

Footplate angle

Set to:

80°	K0040	
90°	K0040	
100°	K0040	

Accessories



1	Tray table <i>(applies to POWER Standing System only)</i>		E0950		\$ 378
	Tray table - hardware		E1028		\$ 316
	Push handles	Pair	K0108		\$ 189
2	Push handles - 2.75" extension	Pair	K0108		\$ 365
3	Push handles - detachable, height adjustable 1" - 6"	Pair	K0108		\$ 378
4	Anti Tips - traditional- Pair	Pair	E0971		\$ 164
5	Tip Aid - attendant- Pair	Pair	K0108		\$ 92
6	Anti Tips - swing away with attendant tip aid	Pair	K0108		\$ 486
	XL Adapter* <i>"Necessary with swing away anti-tips for seat depth positions- of 19.5", 20" & 21" when not using E-Fix or E-Motion wheels- Use of E-Fix or E-Motion wheels require Anti-tips provided by Alber</i>	Pair	K0108		\$ 164
7	Push to Lock- Wheel lock extensions - 2.5"	Pair	K0961		\$ 162
8	Seat depth and footrest quick adjustment kit <i>Includes 8 clamps with post and 4 threaded knobs</i>		K0108		\$ 139
10	Frame reinforcement kit <i>Necessary when e-fix or similar system is mounted Anti-tips - swing away required</i>		K0108		\$ 313
11	Spoke Guards- Transparent	Pair	K0065		\$ 201
	Vehicle tie-down kit		K0108		\$ 1,067

*Fixed attachment points for chair transport Includes frame reinforcement kit
Anti-tips - swing away required- Use of E-Fix or E-Motion wheels require Anti-tips provided by Alber*

Mandatory if client is transported in vehicle. Chair required to be secured with the Q'Straint 5001-T2 or any restraint system certified by Q'Straint 5001-T2. Occupant must be secured using the Q'Straint Vehicle 3 point or any restraint system certified by ISO 10542/2 standards.

Remarks